



# Mesial Distal Sizing Guide

## ANTERIOR REFILLS

<b>E</b>	mm	<b>F</b>	mm	<b>D</b>	mm
E1	6.20	F1	6.20	D1	4.60
E2	6.40	F2	6.40	D2	5.00
E3	6.80	F3	6.80	D3	5.40
E4	7.15	F4	7.15	D4	5.90
E5	7.55	F5	7.55	D5	6.40
E6	7.90	F6	7.90	D6	6.90

  

<b>G</b>	mm	<b>C</b>	mm	<b>H</b>	mm
G1	4.60	C1	6.00	H1	6.00
G2	5.00	C2	6.60	H2	6.60
G3	5.40	C3	7.15	H3	7.15
G4	5.90	C4	7.65	H4	7.65
G5	6.40	C5	8.55	H5	8.55
G6	6.90	C6	9.00	H6	9.00

  

<b>U</b>	mm	<b>SL</b>	mm	<b>SL</b>	mm
U1	4.46	C1 SL	5.36	H1 SL	5.36
U2	4.46	C2 SL	5.90	H2 SL	5.90
U3	4.85	C3 SL	6.39	H3 SL	6.39
U4	4.85	C4 SL	6.83	H4 SL	6.83
		C5 SL	7.64	H5 SL	7.64
		C6 SL	8.04	H6 SL	8.04

## MAXILLARY REFILLS

<b>A</b>	mm	<b>J</b>	mm	<b>B</b>	mm	<b>I</b>	mm
A2	9.15	J2	9.15	B2	6.90	I2	6.90
A3	9.65	J3	9.65	B3	7.30	I3	7.30
A4	10.10	J4	10.10	B4	7.75	I4	7.75
A5	10.40	J5	10.40	B5	8.10	I5	8.10
A6	10.65	J6	10.65	B6	8.50	I6	8.50
A7	11.15	J7	11.15	B7	8.80	I7	8.80
<b>SL</b>	mm	<b>SL</b>	mm	<b>SL</b>	mm	<b>SL</b>	mm
A2 SL	7.91	J2 SL	7.91	B2 SL	6.55	I2 SL	6.55
A3 SL	8.35	J3 SL	8.35	B3 SL	6.94	I3 SL	6.94
A4 SL	8.74	J4 SL	8.74	B4 SL	7.36	I4 SL	7.36
A5 SL	9.00	J5 SL	9.00	B5 SL	7.70	I5 SL	7.70
A6 SL	9.21	J6 SL	9.21	B6 SL	8.00	I6 SL	8.00
A7 SL	9.65	J7 SL	9.65	B7 SL	8.36	I7 SL	8.36

## MANDIBULAR REFILLS

<b>S</b>	mm	<b>L</b>	mm	<b>T</b>	mm	<b>K</b>	mm
S2	7.40	L2	7.40	T2	9.50	K2	9.50
S3	7.80	L3	7.80	T3	10.00	K3	10.00
S4	8.00	L4	8.00	T4	10.40	K4	10.40
S5	8.40	L5	8.40	T5	10.65	K5	10.65
S6	8.90	L6	8.90	T6	11.15	K6	11.15
S7	9.30	L7	9.30	T7	11.60	K7	11.60
<b>SL</b>	mm	<b>SL</b>	mm	<b>SL</b>	mm	<b>SL</b>	mm
S2 SL	6.95	L2 SL	6.95	T2 SL	8.46	K2 SL	8.46
S3 SL	7.33	L3 SL	7.33	T3 SL	8.90	K3 SL	8.90
S4 SL	7.52	L4 SL	7.52	T4 SL	9.26	K4 SL	9.26
S5 SL	7.90	L5 SL	7.90	T5 SL	9.48	K5 SL	9.48
S6 SL	8.40	L6 SL	8.40	T6 SL	9.92	K6 SL	9.92
S7 SL	8.75	L7 SL	8.75	T7 SL	10.32	K7 SL	10.32



**EZCROWNS**

EZCrowns are covered by one or more of US patents Nos. 9,655,690 and D662,594; and multiple other patents pending. Sprig Oral Health Technologies, Inc. © All Rights Reserved 2019.